

AccuSign® BZO

One-Step Benzodiazepines Test

For *In Vitro* Use Only

Simple One-Step Immunoassay for the Qualitative Detection of Benzodiazepines and/or their Metabolites in Human Urine

Catalog No.	DOA-207-35	35 Test Kit
	DOA-207-10	10 Test Kit

Intended Use

AccuSign® BZO is a simple, one-step immunoassay intended for use in the qualitative detection of benzodiazepines in human urine with a cutoff at 300 ng/mL for oxazepam.

The AccuSign® BZO test provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography, mass spectrometry (GC/MS) is the preferred confirmatory method. Other chemical confirmatory methods are available. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.¹

Summary and Explanation

Benzodiazepines are a class of frequently prescribed central nervous system (CNS) depressants which include widely used drugs such as chlordiazepoxide, diazepam, and oxazepam. They have medically useful properties, including antianxiety, sedative, anti-convulsant, and hypnotic effects.² They are taken orally or sometimes by injection, and have a low potential for physical or psychological dependence. Benzodiazepines induce drowsiness and muscle relaxation. Their use can also result in intoxication, similar to drunken behavior except without evidence of alcohol use, and the loss of inhibitions. Chronic abuse can result in addiction and tardive dyskinesia (involuntary muscle movements of the face, limbs, and trunk). Overdose can result in coma and possible death. Withdrawal syndrome includes anxiety, insomnia, tremors, delirium, and convulsions.

The effects of benzodiazepine last 4–8 hours. The different benzodiazepines are absorbed at different rates, and the timing of their psychoactive effects varies with the absorption rate. The drugs are excreted in the urine primarily as the parent compounds or as oxazepam glucuronide, an inactive metabolite (in the case of chlordiazepoxide and diazepam) and are detectable for 1–2 days. Oxazepam is detectable in the urine for up to 7 days.^{2,3}

Principle

The **AccuSign® BZO** test uses solid-phase chromatographic membrane immunoassay technology for the qualitative detection of benzodiazepine. The test is based on the principle of the highly specific immunochemical reactions between antigens and antibodies which are used for the analysis of specific substances in biological fluids. The test relies on the competition between the drug conjugates and the drugs which may be present in the urine sample, for binding to antibodies. In the test procedure, a sample of urine is placed in the Sample well of the device and is allowed to migrate upward. If the drug is present in the urine sample, it competes with the drug conjugate bound to the dye, for the limited antibodies immobilized on the membrane. If the level of drug or drug metabolite is above the cutoff level, the drug will saturate the antibodies, thus inhibiting the binding of the dye coated with drug conjugates to the antibodies on the membrane. This prevents the formation of a line on the membrane. Therefore, a drug-positive urine sample will not generate a line at Test (T) position in the Result window, indicating a positive result from positive drug competition. A negative urine sample will generate a line at Test position in the Result window, indicating a negative result from an absence of competition with free drugs.

In addition to the line that may appear at the Test position in the Result window, a Control line must appear at the Control (C) validation position in the Result window to confirm the viability of the test. This Control line should always be seen if the test is conducted properly. This works as a procedural control, confirming that proper sample volume was used and the reagent system worked. If insufficient sample volume is used, there may not be a Control line, indicating that the test is invalid.


Materials Provided

The **AccuSign® BZO** test kit contains all the reagents necessary to perform the tests.

- **AccuSign® BZO** device. The test device contains a membrane strip coated with polyclonal anti-benzodiazepine antibody and a pad containing drug-dye conjugate in a protein matrix.
- Disposable specimen dispensers.
- Instructions for use.

Precautions

- For *in vitro* diagnostic use only.
- Avoid cross contamination of urine samples by using a new urine specimen container and dropper for each urine sample.
- The test kit does not contain any HIV or hepatitis infective components.
- Urine specimens are potentially infectious. Proper handling and disposal methods should be established according to good laboratory practices.
- The **AccuSign®** device should remain in its original sealed pouch until ready for use. Do not use the test if the pouch is damaged or the seal is broken.
- Do not use the test kit after the expiration date.



Storage and Stability

The **AccuSign® BZO** test kit should be stored at 2–30°C (35–86°F) in the original sealed pouch. The expiration dating was established under these storage conditions.

Specimen Collection and Preparation

Approximately 110 µL of urine sample is required for each test. Fresh urine specimens do not require any special handling or pretreatment. Specimens should be collected in a clean glass or plastic container. If testing will not be performed immediately, specimens should be refrigerated (2–8°C) or frozen. Specimens should be brought to room temperature before testing.

Specimens containing a large amount of particulate matter may give inconsistent test results. Such specimens should be clarified by centrifuging or allowing to settle before testing.

Test Procedure

The test procedure consists of adding the urine sample to the Sample well of the device and watching for the appearance of colored lines in the result window.

Test Protocol

1. For each test, open one **AccuSign® BZO** pouch and label the **AccuSign®** device with the patient ID.
2. Holding the dropper vertically, dispense 3 drops (110 µL) of the urine sample into the Sample well (S).
3. Read the result after 3 minutes, but within 10 minutes of sample addition.

Interpretation of Results

Negative: The appearance of a reddish-purple Control line (C) and a line next to T indicates a negative test result; i.e., no drug above the cutoff level has been detected. The color intensities of the Control line and the Test line may not be equal. *Any faint line in the Result window, visible in 10 minutes, should be interpreted as negative. A negative test result does not indicate the absence of drug in the sample; it only indicates the sample does not contain drug above the cutoff level in qualitative terms.*

Positive: The appearance of only a reddish-purple Control line and no distinct line next to T indicates the test result is positive for BZO (i.e., the specimen contains BZO at a concentration above the cutoff level). *A positive test result does not provide any indication of the level of intoxication or urinary concentration of the drug in the sample; it only indicates the sample contains drug above the cutoff level in qualitative terms.*

Invalid: A distinct Control line (C) should always appear. The test is invalid if no Control line forms at the C position. Such tests should be repeated with a new **AccuSign® BZO** test device.

Limitations

- The test is designed for use with unadulterated urine only.
- There is a possibility that factors such as technical or procedural errors, as well as other substances in the urine sample which are not listed in Table 3 below, may interfere with the test and cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results. If adulteration is suspected, the test should be repeated with a new sample.
- The test result read after 10 minutes may not be consistent with the original reading obtained within the 10 minute reading period. The test must be read within 10 minutes of sample application.

User Quality Control

Internal Control: Each **AccuSign®** test device has built-in controls. The Control line is an internal positive procedural control. A distinct reddish-purple Control line should always appear at the C position, if the test procedure is performed properly, an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the test reagents at the control line and the conjugate-color indicator are reactive. In addition, if the test has been performed correctly and the device is working properly, the background in the result window will become clear and provide a distinct result. This may be considered an internal negative procedural control.

The positive and negative procedural controls contained in each **AccuSign®** test device satisfy the requirements of testing a positive control and a negative control on a daily basis. If the Control line does not appear at the Control position, the test is invalid and a new test should be performed. If the problem persists, contact PBM for technical assistance.

External Control: External controls may also be used to assure that the reagents are working properly and that the assay procedure is followed correctly. It is recommended that a control be tested at regular intervals as good laboratory testing process. For information on how to obtain controls, contact PBM's Technical Services.

Expected Values

AccuSign® BZO is a qualitative assay. The amount of benzodiazepines and/or their metabolites in the urine cannot be estimated by the assay. The assay results distinguish positive from negative samples. Positive results indicate the samples contain benzodiazepine metabolites above the cutoff concentration.

Performance Characteristics

The **AccuSign® BZO** test has been shown to detect oxazepam at an average cutoff of 300 ng/mL in urine. The test also detects other benzodiazepines listed below at the minimum concentrations indicated (Table 2).

The accuracy of **AccuSign® BZO** was evaluated in comparison to a commercially available immunoassay (Syva® EMIT® II) at a cutoff of 300 ng/mL. A total of 223 samples was tested by both procedures. The overall accuracy of the test was 98.7%, as shown below. (Table 1.)

Table 1. Accuracy: Comparison of AccuSign® BZO with Syva® EMIT® II

		Syva® EMIT® II (BZO)		TOTAL
		Positive	Negative	
AccuSign® BZO	Positive	84	2	86
	Negative	1	136	137
TOTAL		85	138	223
		Relative Sensitivity	Relative Specificity	
AccuSign® BZO		98.8% (84/85)	98.6% (136/138)	

Discrepant samples for BZO were analyzed by GC/MS. The one false-negative sample contained the drug at a level of 344 ng/mL, while the two false-positive samples showed 274 and 289 ng/mL.

In a separate study, **AccuSign® BZO** was evaluated against 27 specimens confirmed as positive by GC/MS. The range of drug values was 312 to 610 ng/mL. The results demonstrate the excellent correlation of **AccuSign® BZO** with GC/MS.

Precision and Accuracy

The precision of the **AccuSign® BZO** assay was determined by carrying out the test with serially diluted standard drug solutions using 3 lots of products on 3 different dates. Ninety-eight percent (98%) of the spiked samples containing oxazepam concentrations 25% over the cutoff level (i.e., 375 ng/mL) showed positive results.

Distribution of Random Error:

Twenty (20) blind samples prepared by spiking various concentrations of drug were separately tested by two operators. The test results from the two operators showed complete agreement.

Reproducibility

The reproducibility of the test results of the **AccuSign® BZO** assay was examined at three different sites using a total of 15 blind controls, consisting of 5 negative samples, 5 moderately positive samples (600 ng/mL oxazepam), and 5 strongly positive samples (1,200 ng/mL oxazepam). The results obtained at these three sites with these controls demonstrated 100% agreement with each other.

Specificity

Compounds that are detected by the **AccuSign® BZO** test are listed below (Table 2). The specificity of **AccuSign® BZO** was determined by adding various drugs and drug metabolites to drug-negative urine specimens and testing with the **AccuSign® BZO** test kit. The results are expressed in terms of the concentration required to produce a positive result.

Table 2. Specificity

Compound	Concentration (ng/mL)
Alprazolam	10,000
Bromazepam	2,500
Chlordiazepoxide	1,000
Clobazam	10,000
Clonazepam	10,000
Clorazepate dipotassium	250
Delorazepam	1,000
N-Desalkylflurazepam	1,250
N-Desmethyl diazepam	500
Diazepam	2,000
Estazolam	500
Flunitrazepam	>10,000
7-amino flunitrazepam	4,000
a-Hydroxyalprazolam	250
a-Hydroxytriazolam	2,500
Lorazepam	750
Lormetazepam	8,000
Medazepam	2,500
Midazolam	4,000
Nitrazepam	2,500
Nordiazepam	500
Oxazepam	300
Prazepam	9,000
Temazepam	800
Triazolam	>10,000

The following compounds show no cross-reactivity when tested with **AccuSign® BZO** at a concentration of 100 µg/mL. (Table 3.)

Table 3. Non Cross-Reacting Compounds


Acetaldehyde	ride	amphetamine
4-Acetamidophenol	Ecgonine methyl-ester	Methaqualone
Acetaminophen	(+) Ephedrine	Methoxyphenamine
Acetone	(±) Ephedrine	(±) 3,4-Methylene-dioxyamphet-amine
Acetophenetidin (Phenacetin)	(-) Ephedrine	(±) 3,4-Methylene-dioxymetham-phetamine
N-Acetylprocainamide	(-) Ψ Ephedrine	Methylphenidate
Acetylsalicylic acid	Epinephrine	Methyprylon
Albumin	Erythromycin	Morphine
Aminopyrine	β-Estradiol	Morphine-3-β-D-glucuronide
Amitriptyline	Estriol	Nalidixic acid
Amobarbital	Estrone-3-sulfate	Nalorphine
Amoxapine	Ethyl-p-amino-benzoate	Naloxone
Amoxicillin	Fenopropfen	Naltrexone
D,L-Amphetamine	Furoxime	Naproxen
L-Amphetamine	Gentisic acid	Niacinamide
Ampicillin	Guaiacol glycerol ether	Nifedipine
Apomorphine	Glucose	Norcodein
Ascorbic acid	Glucuronic acid	Norethindrone
Aspartame	Glutethimide	Noroxymorphone
Atropine	Guaiifenesin	D-Norpseudo-ephedrine
Benzilic acid	Hemoglobin	Noscapine
Benzocaine	Hippuric acid	Nylidrin
Benzoic acid	Hydralazine	D,L-Octopamine
Benzoyllecgonine	Hydrochlorothiazide	Oxalic acid
Benzphetamine	Hydrocodone	Oxolinic acid
Bilirubin	Hydrocortisone	Oxycodone
Butabarbital	Hydromorphone	Oxymetazoline
Cannabidiol	O-Hydroxyhippuric acid	Oxymorphone
Cannabinol	3-Hydroxytyramine	Papaverine
Chloralhydrate	Ibuprofen	Penicillin-G
Chloramphenicol	Imipramine	Pentazocaine
Chlorothiazide	Iproniazid	Pentobarbital
Chlorpheniramine	(-) Isoproterenol	Perphenazine
Chlorpromazine	Isoxsuprine	Phencyclidine
Chloroquine	Ketamine	Phendimetrazine
Cholesterol	Ketones	Phenelzine
Clomipramine	Ketoprofen	β-Phenethylamine
Clonidine	Labetalol	Phenobarbital
Cocaine hydrochloride	Levorphanol	Phenothiazine
Codeine	Lidocaine	Phentermine
Cortisone	Loperamide	Phentoin
(-) Cotinine	Loxapine succinate	Phenylbutazone
Creatinine	Lysergic acid diethylamide	L-Phenylephrine
Deoxycorticosterone	Maprotiline	D,L-Phenylpropanol-amine
Dextromethorphan	Melanin	Prednisolone
Dextropropoxyphene	Meperidine	Prednisone
Diclofenac	Meprobamate	Procaine HCl
Diethylpropion	Methadone	
Diflunisal	D-Methamphetamine	
Digoxin	p-Hydroxymeth-	
Diphenhydramine		
Domperidone		
Doxylamine		
Ecgonine hydrochloride		

Promazine	Sodium chloride	D,L-Thyroxine
Promethazine	Sulfamethazine	Tolbutamide
D,L-Propranolol	Sulindac	Triamterene
Propiomazine	Tetracycline	Trifluoperazine
D-Propoxyphene	Tetrahydrocortisone	Trimethoprim
D-Pseudoephedrine	Δ ⁸ -THC	Trimipramine
L-Pseudoephedrine	Δ ⁹ -THC	Tryptamine
Quinidine	11-nor-Δ ⁹ -THC-9-COOH	D,L-Tryptophan
Quinine	COOH	Tyramine
Rantidine	Tetrahydrozoline	D,L-Tyrosine
Salicylic acid	Thebaine	Uric acid
Secobarbital	Thiamine	Verapamil
Serotonin	Thioridazine	Zomepirac

References

- Hawks RL, Chiang CN, eds. *Urine Testing for Drugs of Abuse*. Rockville, MD: National Institute on Drug Abuse (NIDA), Research Monograph 73;1986.
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 2nd Ed., Davis, CA: Biomedical Publ.;1982.
- Greenblatt DJ, Shader RI. *Benzodiazepines in Clinical Practice*. New York: Raven Press;1974.

Symbols Key

	Manufactured by
	CE Mark
	Authorized Representative
	In Vitro Diagnostic Medical Device
	Catalog Number
	Consult Instructions for Use
	Batch Code
	“Use By” date in year-month-day format
	Temperature Limitation
	Contains sufficient for <n> tests
	Do not reuse
	Contents
	Test Device
	Transfer Pipette
	Instructions for Use
	One-step immunochromatographic Assay for the Detection of Drugs of Abuse in Urine
	Benzodiazepines Test

AccuSign® is a Registered Trademark of Princeton BioMeditech Corporation.

Patent No.: 5,559,041



© 2000 PBM
Printed in U.S.A.
Revised Oct 2003
P-5846-D 1001BL



MT Promedt Consulting GmbH
Eisenbahnstrasse 2
D-66386 St. Ingbert
Germany
+49-68 94-58 10 20



Princeton BioMeditech Corporation
Princeton, NJ 08543-7139 U.S.A.
1-732-274-1000 www.pbmc.com