

Non Cross-Reacting Compounds

4-Acetaminophenol
Acetophenetidin
Albumin
Amitypyline
Amoxapine
Ampicillin
Aspartame
l-Amphetamine
Benzoic acid
Benzoylcocgonine
B-Estradiol
Brompheniramine
(+)-Chlorpheniramine
Caffeine
Cannabinol
Chloramphenicol
Chloroquine
Chlorprothixene
Cimetidine
Clonidine
Codeine
Creatinine
Cyclobenzaprine
4-Dimethylaminoantipyrine
Deoxycorticosterone
Diazepam
Digoxin
Disopyramide
R(-)Deprenyl
(±)Epinephrine
Ecgonine
EDDP
EMDP
Erythromycin
Ethanol
Etodolac
l-Epinephrine
Fenfluramine
Fentanyl
Furosemide
Gentisic acid
3-Hydroxytyramine
Hydralazine
Hydrocodone
Hydromorphone
o-Hydroxyhippuric acid
p-Hydroxyamphetamine
(-)isoproterenol
Imipramine
Isosuprine
Ketamine
Labetalol
Lidocaine
Lithium Carbonate
d-Methamphetamine
Maprotiline
MDMA**
Mepentermine
Methadone
Methoxyphenamine
Metoprolol
Morphine Sulfate
a-Naphthaleneacetic acid
d-Norpropoxyphene
Nalorphine
Naltrexone
Norcodeine
Normorphine
D,l-Octopamine
Oxalic acid
Oxilinic acid
Oxymetazoline
(±)Phenylpropanolamine
D,l-Propranolol
d-Pseudoephedrine
Papaverine
Penicillin-G
Pentobarbital
Phencyclidine
Pheniramine

Acetone
Acetylsalicylic Acid
Aminopyrine
Amobarbital
Amoxicillin
Apomorphine
Atropine
l-Ascorbic acid
Benzoic acid
Benzphetamine
Bilirubin
Buspirone
(±)Chlorpheniramine
Cannabidiol
Chloral Hydrate
Chlordiazepoxide
Chlorpromazine
Cholestrol
Clomopramme
Cocaine
Cortisone
Cyclobarbitol
(-)Deoxyephedrine
5,5Diphenylhydantoin
Dextromethorphan
Diclofenax
Diphenhydramine
Doxylamine
(-)ψ-Ephedrine
(tr,2s)-(-)Ephedrine
Ecgonine Methyl ester
Efavirenz (Sustiva)
Ephedrine
Estrone 3-sulfate
Ethyl-p-aminobenzoate
l-Ephedrine
Famprofazone
Fenoprofen
Fluoxetine
d(-)Glucose
Guaiaicol Glyceryl Ether
Hemoglobin
Hydrochlorothiazide
Hydrocortisone
Hydroxyzine
p-Hydrizymethamphetamine
p-Hydroxynorephedrine
l-bupropfen
lproniazid
Kanamycin
Ketoprofen
Levorphanol
Lindane
Loperamine
l-Methamphetamine
MDA*
Meperidine
Meprobamate
Methaqualone
Methyprylon
Morphine 3-β-d-glucuronide
Mythlphenidate
d-l Norephedrine
Nalidixic acid
Naloxone
Nimesulide
Norethindone
Noscapine
Orphenadrine
Oxazepam
Oxycodone
Oxymorphone
B-Phenylethylamine
d-Propoxyphene
l-Phenylephrine
Pemoline
Pentazocine
Perphenazine
Phenelzine
Phenobarbital

Phenothiazine
Prednisolone
Procaine
Protrazine
Quinacrine
Quinine
Riboflavin
Secobarbital
Sodium Chloride
Sulindac
Tetracycline
3-acetate Tetrahydrozoline
Theophylline
Buspirone
Thiodiazine
Trifluoperazine
Thrumethoprim
d,l-Tryptophan
D,l-Tyrosine
Verapmil

Phentermine
Prednisone
Promethazine
Trans-2-phenylcyclopropylamine
Quinidine
Ramtidine
Salicylic acid
Serotonin
Sulfamethazine
Temazepam
Tetrahydrocortisone
Thebaine
Thamine
l-Thyroxine
Trimethobenzamide
Triprimamine
Tyramine
Uric Acir
Zomepirax

*MDA = 3,4-Methylenedioxyamphetamine
**MDMA = 3,4-Methylenedioxyamphetamine

REFERENCES

1. U.S. Department of Health and Human Services, Public Health Services, Centers for Disease Control: *Reducing the Health Consequences of Smoking, 25 Years of Progress, A Report of the Surgeon General*, Rockville, MD. Office on Smoking and Health, 1989
2. U.S. Department of Health and Human Services, Public Health Service, *The Health Consequences of Smoking: Cardiovascular Disease, A Report of the Surgeon General*, Rockville, MD. Office on Smoking and Health 1983
3. Bruckert E, Jacob N, Lamaire L, Truffert J, PercheronF, and de Gennes JL. *Relationship between Smoking Status and Serum Lipid in a Hyperlipidemic Population and Analysis of Possible Confounding Factors*. Clin Chem 1992;38:1698-1705.
4. Projer. R, Whitfield JB, Poulos V, Eckhard IF, Richmond R, Hensley WJ. *Carboxyhaemoglobin, Cotinine and Thiocyanate Assay Compared for Distinguishing Smokers from Non-smokers*. Clin Chem 1984;30:1377-1380
5. Benowitz NL, Kuyt F, Jacob P, et al. *Cotinine Disposition and effects*. Clin Pharmacol Ther 34, 604-611 (1983).
6. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 6th Edition Biomedical Publications, Foster City, CA. 2002; 744-747
7. Hardman JG and Limbird LE, Goodman and Gilman's: *The Pharmacological Basis for Therapeutics*. 10th Edition. McGraw Hill Medical Publishing, 2001; 208-209

DTA Pty Ltd, South Africa
PO Box 1661, Sun Valley 7985 Cape Town
email@drugtesting.co.za
www.drugtesting.co.za



AccuTest™ NIC Test Cassette

FOR THE QUALITATIVE ASSESSMENT OF NICOTINE
AND THE
METABOLITES IN HUMAN URINE



For in vitro Diagnostic and Forensic Use

INTENDED USE

The AccuTest™ NIC test device is an immunochromatography based one step in vitro test. It is designed for qualitative determination of drug substances in human urine specimens. This assay may be used in the point of care setting. The cut-off concentration using our test is as follows;

Nicotine 200 µg/ml of cotinine

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. The Substance Abuse Mental Health Services Administration (SAMHSA) has established gas chromatography/ mass spectrometry (GC/MS) as the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

Note: It is not for use in diagnosing disease or illness; it is for investigational use only.

SUMMARY AND EXPLANATION

Smoking has been identified as a major risk factor for lung cancer, and cardiovascular disease.^{1,2} Self-reporting of smoking status is not reliable.³ The determination of cotinine, a major metabolite of nicotine, has become the preferred biomedical method of assessing the smoking status of individuals on account of its sensitivity and specificity.⁴ Cotinine is present in blood, urine and saliva of individuals who smoke or chew tobacco or who inhale tobacco smoke produced by others. As an objective indicator of nicotine intake or confirmation of nonsmoker status, cotinine offers several advantages over other biochemical measures: it is a specific indicator of nicotine intake, its concentrations are not influenced by confounding factors such as diet or environment, its average biological half-life in blood is 19 hours, and its concentration within a given individual varies by on 15 to 20% over the course of a day.⁵ Cotinine assay is thus a superior objective measure of exposure to nicotine.

PRINCIPLE

The AccuTest™ NIC test device is based on the principle of specific immunochemical reaction between antibodies and antigen to analyze particular compound in human urine specimen. The assay relies on the competition for binding antibody. When drug is present in the urine specimen, it competes with drug conjugate for the limited amount of antibody dye conjugate. When the amount of drug is equal or more than the cut-off, it will prevent the binding of drug conjugate to the antibody. Therefore, a positive urine specimen will not show a colored band on the test line zone, indicating a positive result, while the presence of a colored band indicates a negative result. A control line is present in the test window to work as procedural control. This colored band should always appear on the control line zone if the test device is stored in good condition and the test is performed appropriately.

MATERIAL PROVIDED

1. An AccuTest™ NIC test device (cassette).
The amount of each coated antigen and/or antibody on the strip is less than 1.0 mg for antigen conjugate and is less than 1.0 mg for goat anti-mouse IgG antibody.
Test zone: contains drug bovine protein antigen conjugates

- Control zone: contains Goat anti-mouse IgG antibody
 Conjugate pad: contains mice monoclonal anti-drug antibody.

2. Dropper
3. Instruction for use.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Urine collection container.
2. Timer or clock.

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic and professional use only.
2. Do not use the test device beyond the expiration date.
3. Urine specimens may be infectious; properly handle and dispose of all used reaction devices in a biohazard container.
4. Visually inspect the foil package to insure it is intact. If the package is not intact, the integrity of the device might be compromised.
5. Use a new urine specimen cup for each sample to avoid cross contamination.

STORAGE AND STABILITY

The test device should be stored at 2°C to 28°C; **do not freeze** and will be effective until the expiration date stated on the package. The product is humidity-sensitive and should be used immediately after being open. Any improperly sealed product should be discarded.

SPECIMEN COLLECTION AND PREPARATION

Fresh urine does not require any special handling or pretreatment. Specimen should be collected in a clean, dry, plastic or glass container. If the assay is not performed immediately, urine specimen may be refrigerated at 2- 8 °C or frozen up to 7 days. Specimens should be brought to room temperature before testing. Urine specimens exhibiting a large amount of precipitate or turbidity should be centrifuged or allowed to settle before testing. Avoid contact with skin by wearing gloves and proper laboratory attire.

QUALITY CONTROL

Good Laboratory practice recommends the daily use of control materials to validate the reliability of device. Control materials should be assayed as clinical specimen and challenging to the assay cutoff concentration, e.g., 25% above and below cutoff concentration. If control values do not fall within establish range, assay results are invalid. Control materials that are not provided with this test kit are commercially available.

The AccuTest™ NIC urine screen test device provides a built-in process control with a different antigen/antibody reaction at the control region (C). This control line should always appear regardless the presence of drug or metabolite. If the control line does not appear, the test device should be discarded and the obtained result is invalid. The presence of this control band in the control region serve as 1) verification that sufficient volume is added, 2) that proper flow is obtained.

PROCEDURE

1. Bring all materials and specimens to room temperature.
2. Remove the test strip from sealed foil pouch.
3. Place the sample pad end of the test strip (up to the bottom blue line mark) into the urine specimen.
4. Hold the strip in the urine for 10 seconds, remove from the urine.
5. At 5 minutes, read the test strip.

Do not interpret the result after 5 minutes. Waiting more than five minutes may cause the reading to be inaccurate. To avoid confusion, discard the test strip after interpreting the result.

INTERPRETATION OF RESULTS

Negative:

Two colored bands form on any strip of the card. The appearance of two colored bands, one in test line zone and the other in control line zone, indicates negative result for that particular test(s). The negative result does not indicate the

absence of drug in the specimen; it only indicates the level of tested drug in the specimen is less than cut-off level.


Positive:

One colored band appears in control line zone. No coloured band is found in test line zone. This is an indication the level of tested drug(s) in the specimen is above the cut- off level.

Invalid:

If there is no colored band in control line zone of any strip, the test result is invalid. Retest the sample with a new device.

Note: A borderline(+/-) in test line zone should be considered negative result.



LIMITATION OF PROCEDURE

The assay is designed for use with human urine only. A positive result with any of the tests indicates only the presence of a drug/metabolite and does not indicate or measure intoxication. There is a possibility that technical or processed urine error as well other substances in certain foods and medicines may interfere with the test and cause false results. Please refer "SPECIFICITY" section for lists of substances that will produce either positive results, or that do not interfere with test performance. If a drug/metabolite is found present in the urine specimen, the assay does not indicate frequency of drug use or distinguish between drug of abuse and certain foods and medicines.

EXPECTED RESULTS

The AccuTest™ NIC test device is a qualitative assay. It identifies the drug in human urine at its cut-off concentration or higher. The concentration of the drug cannot be determined by this assay. The test is intended to distinguish negative result from presumptive positive result. All positive results must be confirmed using an alternate method, preferably GC/MS.

PERFORMANCE CHARACTERISTICS

Accuracy:

A side-by-side comparison was conducted using the AccuTest™ NIC test device (urine) and a leading commercially available COT rapid test. Testing was performed on 300 clinical specimens collected from smoking and non-smoking volunteers. The following results were tabled;-

METHOD	Other COT Device		TOTAL RESULTS
	Positive	Negative	
AccuTest™ NIC Test	Results		
	Positive	103	12
	Negative	0	185
Total Results		103	197
% Agreement		>99%	94%

Analytical Sensitivity

A drug-free urine pool was spiked with Cotinine at the following concentrations: 0µg/mL, 100 µg/mL, 200 µg/mL, 250 µg/mL, 300 µg/mL and 400 µg/mL. The result demonstrates >99% accuracy at 100% above and 50% below the cut-off concentration. The data is summarized within the following table:-

Cotinine Concentration (µg/mL)	Percent of Cut-off	n	Visual Results	
			Negative	Positive
0	0	90	90	0
100	-50%	90	90	0
150	-25%	90	90	0
200	Cut-off	90	63	27
250	+25%	90	40	50
300	+50%	90	16	74
400	+100%	90	0	90

Analytical Sensitivity

The following table lists compounds that are positively detected in urine by the AccuTest™ NIC test device (urine) at 5 minutes

Compound	Concentration (µg/mL)
(-)-Cotinine	200
(-)-Nicotine	6,500

Precision

A study conducted by trained operators using 2 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing no Cotinine, 50% Cotinine below cut-off level and 100% Cotinine above the 200µg/mL cut-off level were used. The following results were tabled:-

Cotinine Concentration (µg/mL)	n per lot	LOT A		LOT B	
		-	+	-	+
0	30	30	0	30	0
100	30	30	0	30	0
400	30	0	30	0	30

Effects of Urinary Specific Gravity

Fifteen urine specimens of normal, high and low specific gravity ranges were spiked with 100 µg/mL and 400 µg/mL. The AccuTest™ NIC test device (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urine specific gravity do not affect the test results.

Effects of Urinary pH

The pH of an aliquot negative urines pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Cotinine to 100 µg/mL and 400 µg/mL. The spiked, pH-adjusted urine was tested with the AccuTest™ NIC test device (urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test device.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Cotinine positive urine. The following compounds show no cross-reactivity when tested with the AccuTest™ NIC test device at a concentration of 100µg/mL.