


AccuTest™

International Gold Standard for drug screening

Multi Panel Drug Screening Device in Cassette

FOR THE QUALITATIVE ASSESSMENT OF DRUGS AND THEIR METABOLITES
IN HUMAN URINE



- Easy to use
- Accurate and reliable
- Instant results



- FDA APPROVED
- CLIA WAIVED
- CE MARKED

All batches independently tested to assure the products meet our on-going high quality standards



The preferred supplier of high quality sophisticated and innovative diagnostic assays

INTENDED USE


The multi-panel drug screening device in a cassette is a one-step in vitro test device consisting of drugs of abuse screen test strips. It is designed for qualitative determination of drug substances in human urine specimens. This assay may be used in the point of care setting.




Remove cap



Dip into urine sample for 10 seconds



Replace cap



Wait 5 minutes to interpret results

Available from 2 up to 12 drug analyses combinations

Drug Analyses Choices:

Parameter	Callibrator	Cut-off (ng/mL)										
		10	25	50	100	150	200	300	500	600	1000	2000
AMP (Amphetamine)	Amphetamine							✓	✓			✓
BAR (Barbiturates)	Secobarbital							✓				
BZO (Benzodiazepines)	Oxazepam				✓		✓	✓				
BUP (Buprenorphine)	Buprenorphine-β-3-D-glucuronide	✓										
COC (Cocaine)	Benzoylgonine				✓		✓	✓				
COT (Nicotine)	Cotinine						✓			✓		
EDDP (Methadone metabolite)	2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine				✓							
FYL (Fentanyl)	Fentanyl	✓										
KET (Ketamine)	Ketamine											✓
MDA (Methylenedioxyamphetamine)	Methylenedioxyamphetamine								✓			
MDMA (Ecstasy)	3,4-Methylenedioxymethamphetamine							✓	✓			✓
MTD (Methadone)	Methadone							✓				
MET (Methamphetamine)	Methamphetamine							✓	✓			✓
MOR (OPI) (Morphine/Heroin)	Morphine				✓		✓	✓				✓
MQL (Mandrax)	Methaqualone							✓				
OXY (Oxycodone)	Oxycodone				✓							
PCP (Phencyclidine)	Phencyclidine		✓									
PPX (Propoxyphene)	Norpropoxyphene							✓				
TCA (Tricyclic antidepressants)	Nortriptyline											✓
THC (Marijuana)	11-nor-Δ ⁹ -THC-9-COOH		✓	✓		✓	✓		✓			
TML (Tramadol)	Tramadol				✓							



Drug Testing Africa

+27 21 785 5848

+27 86 518 2210

info@drugtesting.co.za

www.drugtesting.co.za