

# AccuTest™

International Gold Standard for drug screening

## Multi Panel Drug Screening Device in Split Sample Collection Container

FOR THE QUALITATIVE ASSESSMENT OF DRUGS AND THEIR METABOLITES AND THE  
LEVELS OF pH AND CREATININE IN HUMAN URINE



- FDA APPROVED
- CLIA WAIVED
- CE MARKED



- Easy to use
- Accurate and reliable
- Instant results


All batches independently tested to assure the products meet our on-going high quality standards



The preferred supplier of high quality sophisticated and innovative diagnostic assays

## INTENDED USE

The multi-panel drug screening device in split sample collection container is a one-step in vitro test device consisting of a drugs of abuse screen test and adulteration test strips. It is designed for qualitative determination of drug substances in human urine specimens. A temperature strip affixed to the container measures the temperature insuring the specimen has not been substituted. Adulteration test strips (pH and creatinine) is built into the test panel to validate the integrity of the urine specimen examined. This assay may be used in the point of care setting.



### Drug Analyses Choices:

Parameter	Callibrator	Cut-off (ng/mL)										
		10	25	50	100	150	200	300	500	600	1000	2000
AMP (Amphetamine)	Amphetamine							✓	✓		✓	
BAR (Barbiturates)	Secobarbital							✓				
BZO (Benzodiazepines)	Oxazepam				✓		✓	✓				
BUP (Buprenorphine)	Buprenorphine-β-3-D-glucuronide	✓										
COC (Cocaine)	Benzoylcegonine				✓		✓	✓				
COT (Nicotine)	Cotinine						✓			✓		
EDDP (Methadone metabolite)	2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine				✓							
FYL (Fentanyl)	Fentanyl	✓										
KET (Ketamine)	Ketamine										✓	
MDA (Methylenedioxyamphetamine)	Methylenedioxyamphetamine								✓			
MDMA (Ecstasy)	3,4-Methylenedioxymethamphetamine							✓	✓		✓	
MTD (Methadone)	Methadone							✓				
MET (Methamphetamine)	Methamphetamine							✓	✓		✓	
MOR (OPI) (Morphine/Heroin)	Morphine				✓		✓	✓				✓
MQL (Mandrax)	Methaqualone							✓				
OXY (Oxycodone)	Oxycodone				✓							
PCP (Phencyclidine)	Phencyclidine		✓									
PPX (Propoxyphene)	Norpropoxyphene							✓				
TCA (Tricyclic antidepressants)	Nortriptyline										✓	
THC (Marijuana)	11-nor-Δ <sup>9</sup> -THC-9-COOH		✓	✓		✓	✓		✓			
TML (Tramadol)	Tramadol				✓							

## Drug Testing Africa

+27 21 785 5848

+27 86 518 2210

[info@drugtesting.co.za](mailto:info@drugtesting.co.za)

[www.drugtesting.co.za](http://www.drugtesting.co.za)

